

DNA Extraction Protocol from 25mm 0.2um filters

Low-throughput, using DNeasy columns

by V. Rich 2/10/07, used for extracting DNA for 454 sequencing

This protocol is adapted from the 96-well protocol I worked out for environmental samples, which is built on the lab Steripak filter extraction protocol and the Suzuki et al 2001 extraction protocol.

Lysis Buffer: (make fresh because of sucrose)
adapted from Steripak protocol, with addition of RNase

<u>Final Concentration</u>	<u>For 20ml</u>
40mM EDTA	1.6ml of 0.5M EDTA
50mM Tris (pH 8.3)	1.0ml of 1M Tris (pH 8.3)
0.73M Sucrose	5.13g of Sucrose

filter sterilize through 0.2um filter and split into 2 10ml aliquots:

1. With Lysozyme & RNase: right before use, add to one 10ml aliquot:

<u>Final Concentration</u>	<u>For 10ml</u>
1.15mg/ml Lysozyme	11.54ug lysozyme
200ug/ml RNase 100mg/ml	20ul RNase 100mg/ml

shake to dissolve thoroughly, then filter-sterilize again

2. With ProK:

weigh out minimum amount ProK, then add the appropriate amount of lysis buffer (from the second 10ml aliquot)

<u>Final Concentration</u>	<u>For e.g. 10ml</u>	<u>For e.g. 600ul</u>
10mg/ml	100mg ProK	6mg ProK

filter sterilize

Protocol

Increasing the lysis volume to try to maximize extraction efficiency, and because I'm using the DNeasy columns not the 96-well plate, so can load larger volumes...

Step 1: Cell Lysis & RNA removal

- thaw filters on ice
- transfer each filter to screw-top, O-ringed eppendorf tube
- add 500ul lysis buffer to each tube
- **incubate 37°C for 30"**, rotating end-over-end at angle, for optimal mixing with minimal frothing

Step 2: Protein Degradation

adapted from lab Steripak protocol

- add 37.5ul of Proteinase K solution (10mg/ml made up in lysis buffer) to a final conc. of 0.65 mg/ml
- add 59.7ul 10% SDS to a final conc of 1%
- **incubate at 55°C for 2 hours**, rotating end-over-end at angle

Step 3: DNA Purification Through Qiagen DNeasy Tissue Kit DNA-binding column

- add 600ul Buffer AL (=Buffer AL/E without the ethanol added)
mix thoroughly by vortexing
incubate 70°C for 10"
- add 600ul 96-100% EtOH
mix by vortexing vigorously
check pH of lysate, must be pH <7 to get max binding efficiency
- pipet max possible onto spin columns
spin >6000 x g for 10"; discard flow-through
- pipet additional lysate aliquots onto spin columns
spin >6000 x g for 1"; discard flow-through and collection tube, transfer column to new collection tube
- Add 500ul Buffer AW1
spin >6000 x g 1"; discard flow-through
- Add 500ul Buffer AW2
spin >6000 x g for 3"; discard flow-through and collection tube, transfer column to new collection tube
- then, to dry columns, spin 20,000x g for 3"
- add 200ul pre-heated 70°C Buffer AE or water
incubate 1" @ RT
spin >6000 x g 2" to elute
- repeat with a second 200ul, this will increase yield up to 25%

Step 4: Final DNA Clean-up by Size-Exclusion Columns

- transfer the eluted DNA to the Microcon 30
- rinse DNA with 500ul Ambion water
- add 20ul dilute TE, pipette up and down 20 times and transfer to a clean 96-well plate for storage
(optional: repeat with another 20ul to ensure all retrieved)